Note: This Class B digital apparatus complies with Canadian ICES-003. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation of the device.

Five (5) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding batteries) for five (5) years for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. For product support & warranty information: Please contact us with proof of purchase (purchase receipt) at:

**Springfield Instruments Inc./
Les Instruments Springfield
104 Barr Street
St.-Laurent, Quebec
H4T 1Y4
1-844-687-4040**

There are no express warranties except as listed above.

©2017 Taylor Precision Products. Taylor® is a registered trademark of Taylor Precision Products. All rights reserved.

**Made to our exact specifications in China.
Not legal for trade.**

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Model # 7422EF

Electronic Scale

Thank you for purchasing a Taylor® electronic scale. All of our precision instruments for weight measurement are engineered and designed to meet the highest quality standards ... to assure you years of uncompromising accuracy and consistently dependable, convenient performance.

Please read this instruction manual carefully before use. Keep these instructions handy for future reference.

General Operation

1. This scale operates on 3 "AAA" batteries, which are included. Please see "Battery Installation" section of this manual for installation instructions. Also, some models have a static cling label on the lens to prevent scratching. Please remove before use.
2. Your scale has been set at the factory to read weight in pounds (lb). You may select kilograms (kg) or stone (st lb) (1 stone = 14 pounds) units of measurement by pressing the "kg-lb-st" button on the bottom of the scale.
3. Place scale on flat, hard surface. Carpeted or uneven floors may affect accuracy.
4. Step on the scale (disregarding the reading) and step off. This process insures the scale is set and is now ready for use. Do this whenever the scale is moved. At all other times step straight onto the scale.
5. Step onto the scale platform and remain still. The display will count up to your weight. The display will flash twice when the weight is locked. Your weight is now displayed.
6. The scale will turn off automatically.

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Battery Installation

- Use 3 AAA batteries (included. Remove any plastic wrap from the batteries before proceeding). Always replace all batteries at the same time; do not combine old and new batteries. Do not mix Alkaline, carbon-zinc (standard) or Nickel-Cadmium (rechargeable) batteries.
- Open the battery compartment on the bottom of the scale.
- Place the new batteries into the battery compartment according to the polarity illustration in the compartment.
- Replace the battery compartment cover.

NOTE: Please recycle or dispose of batteries per local regulations.

WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.

PRECAUTION: Do not dispose of batteries in fire. Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period of time.

Low Battery

Your Taylor® scale is equipped with a low battery indicator. It will display “Lo” or readings will grow dim or become irregular if the batteries need replacement.

Trouble Shooting

In Case of Difficulty...

- Place your feet in the same position every time you weigh yourself. Stand still during measurement.
- Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
- For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
- The maximum capacity of this scale is 440 lbs / 200 kgs / 31 st 6 lbs. Should “Err ” appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.

Precautions

- Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
- Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
- This scale is designed to read weights between 11 and 440 lbs / 5 and 200 kgs / 11 and 31 st 6 lbs. To prevent damage, do not place more than 440 lbs / 200 kgs / 31 st 6 lbs on it.
- Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
- This scale is an extremely sensitive weighing device. In order to prevent run down of the batteries, do not store anything on the scale.
- Do not dispose of batteries in fire. Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period of time.

Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals.

Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

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