Electronic Scale

Leading the Way in Accuracy® Instruction Manual

Five (5) Year Limited Warranty

This scale is warranted against defects in materials of workmanship (excluding batteries) for five (5) years for the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair. Do not return to retailer. For product support & warranty information: Please contact us with proof of purchase (purchase receipt) at:

Springfield Instruments Inc./ Les Instruments Springfield
104 Barr Street
St.-Laurent, Quebec
H4T 1Y4
1-844-687-4040
info@springfieldinstruments.com

There are no express warranties except as listed above.

©2017 Taylor Precision Products Taylor® is a registered trademark of Taylor Precision Products, Inc. All rights reserved.

Made to our exact specifications in China.
Not legal for trade.

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
Battery Replacement

1. Use 1- 3 volt CR2032 lithium battery (as installed).
2. Remove screw (or screws) from the battery compartment cover on the bottom of the scale and remove the cover. Take care when removing and replacing the battery door screw. Do not turn the screw in the wrong direction or over tighten it as this may strip the screw threads in the plastic battery door.
3. Remove the old battery from the compartment.
4. Place the new battery with the "+" sign up into the battery compartment.
5. Replace the battery compartment cover and screw(s).

NOTE: Please recycle or dispose of batteries per local regulations.

WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.

PRECAUTION: Do not dispose of batteries in fire. Batteries may explode or leak. Remove the batteries if the scale will not be used for a long period of time.

Low Battery

Your Taylor® scale is equipped with a low battery indicator. It will display "LO" or readings will grow dim or become irregular if the battery needs replacement.

Trouble Shooting

In Case of Difficulty...

1. Place your feet in the same position every time you weigh yourself. Stand still during measurement.
2. Weigh yourself without clothing whenever possible. Clothes can vary in weight.
3. Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
4. For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
5. The maximum capacity of this scale is 330 lbs / 150 kgs / 24 st. Should "Err" appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.
6. "Err" also indicates an initialization or instability error. Step off the scale. Ensure the scale is on a flat, hard surface. Press on the platform. The scale will turn on and off. Repeat measurement.

Precautions

1. Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.
2. Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
3. This scale is designed to read weights between 11 and 330 lbs / 5 and 150kgs / 2 and 24 st. To prevent damage, do not place more than 330lbs / 150 kgs /24 st on it.
4. Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
5. This scale is an extremely sensitive weighing device. In order to prevent run down of the battery, do not store anything on the scale.
6. Do not dispose of batteries in fire. Batteries may explode or leak. Remove the battery if the scale will not be used for a long period of time.

Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals. Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

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